

FEATURES



Download the HP Hearing App (from the App Store or Google Play) to wirelessly control your hearing aids**



Ear ID™ measures your hearing levels and delivers personalized hearing amplification calibrated to your unique hearing profile



Streams music and audio with high fidelity from Bluetooth® connected devices



Control the noise around you with Active Noise Cancellation



Focuses on sounds in front of you using directional microphones



Each hearing aid is powered by a rechargeable battery for up to 8 hours of battery life and included charge case provides added on-the-go power



Water and dust resistant IP54 rating

2-year manufacturer's repair warranty

www.hphearingpro.com



HP Hearing PRO
OTC Hearing Aids



FDA CLEARED
Medical Device

30% Better Speech Understanding in Noise*
Self-Fitting OTC Hearing Aids with App
Active Noise Cancellation
Bluetooth® Streaming



HP is a trademark of HP Inc. or its affiliates, used by Nuheara under license from HP Inc. Ear ID™ is a trademark of Nuheara Limited registered in Australia, the USA, and other countries.

Google Play and the Google Play logo are trademarks of Google LLC.

Apple, the Apple logo, iPhone, and iPad are trademarks of Apple Inc., registered in the U.S. and other countries and regions. App Store is a service mark of Apple Inc.

Qualcomm aptX is a product of Qualcomm Technologies, Inc. and/or its subsidiaries. Qualcomm is a trademark of Qualcomm Incorporated, registered in the United States and other countries. aptX is a trademark of Qualcomm Technologies International, Ltd., registered in the United States and other countries.


©2022, Hearing Components. Comply is a registered trademark of Hearing Components, Inc.

This information and other labelling, including the user instructional brochure, are available on the internet at: www.hphearingpro.com/labelingandresources

You may also call (844) 540-0384 or write to info@hphearingpro.com or 13555 S.E. 36th St. Suite 100, Bellevue, Washington, USA 98006 to request a paper copy of this information and other labelling.

*Nuheara 510(k) clearance, K221064, for self-fitting, OTC hearing aid.

**Smartphone required which operates iOS 11 or higher or Android which operates Android 10 or higher.



Nuheara Limited
190 Aberdeen Street Northbridge
WA 6003 Australia
Made in Malaysia

Sleeve Part Number PKG0000057-02 Rev 3.1



HP is a trademark of HP Inc. or its affiliates, used by Nuheara under license from HP Inc.

Ear ID™ is a trademark of Nuheara Limited registered in Australia, the USA, and other countries.

Google Play and the Google Play logo are trademarks of Google LLC.

Apple, the Apple logo, iPhone, and iPad are trademarks of Apple Inc., registered in the U.S. and other countries and regions. App Store is a service mark of Apple Inc.

Qualcomm aptX is a product of Qualcomm Technologies, Inc. and/or its subsidiaries. Qualcomm is a trademark of Qualcomm Incorporated, registered in the United States and other countries. aptX is a trademark of Qualcomm Technologies International, Ltd., registered in the United States and other countries.

© 2022, Hearing Components. Comply is a registered trademark of Hearing Components, Inc.

This information and other labelling, including the user instructional brochure, are available on the internet at:

www.hphearingpro.com/labelingandresources

You may also call (844) 540-0384 or write to info@hphearingpro.com or 13555 S.E. 36th St. Suite 100, Bellevue, Washington, USA 98006 to request a paper copy of this information and other labelling.

*Nuheara 510(k) clearance, K221064, for self-fitting, OTC hearing aid.

**Smartphone required which operates iOS 11 or higher or Android which operates Android 10 or higher.




Nuheara Limited
190 Aberdeen Street Northbridge
WA 6003 Australia
Made in Malaysia

Sleeve Part Number PKG0000057-02 Rev 3.1



In The Box



This hearing aid is for adults with signs of mild to moderate hearing loss. How do you know if you have this?

- You have trouble hearing speech in noisy places
- You find it hard to follow speech in groups
- You have trouble hearing on the phone
- Listening makes you tired
- You need to turn up the volume on the TV or radio, and other people complain it's too loud

Some people with hearing loss may need help from a hearing healthcare professional. How do you know if you need to see one?

- You can't hear speech even if the room is quiet
- You don't hear loud sounds well, for example, you don't hear loud music, power tools, engines, or other very noisy things

If your hearing loss makes it hard to hear loud noises, this hearing aid may not be your best choice without help from a professional. If this hearing aid does not help you enough, ask for help from a hearing healthcare professional.

WARNING: When to See a Doctor
If you have any of the problems listed below, please see a doctor, preferably an ear-nose-throat doctor (an ENT).

- Your ear has a birth defect or an unusual shape. Your ear was injured or deformed in an accident
- You saw blood, pus, or fluid coming out of your ear in the past 6 months
- Your ear feels painful or uncomfortable
- You have a lot of ear wax, or you think something could be in your ear
- You get really dizzy or have a feeling of spinning or swaying (called vertigo)
- Your hearing changed suddenly in the past 6 months
- Your hearing changes: it gets worse then gets better again
- You have worse hearing in one ear
- You hear ringing or buzzing in only one ear

WARNING: If you are younger than 18, do not use this.

You should go to a doctor, preferably an ear-nose-throat doctor (an ENT), because your condition needs specialized care. Over-the-counter hearing aids are only for users who are age 18 or older.

If you purchased your product anywhere other than www.hphearingpro.com, then you must return it to the original place of purchase for a refund or exchange, observing the original place of purchase return, exchange, and money-back policies.

FEATURES



Download the HP Hearing App (from the App Store or Google Play) to wirelessly control your hearing aids**



Ear ID™ measures your hearing levels and delivers personalized hearing amplification calibrated to your unique hearing profile



Streams music and audio with high fidelity from Bluetooth® connected devices



Control the noise around you with Active Noise Cancellation



Focuses on sounds in front of you using directional microphones



Each hearing aid is powered by a rechargeable battery for up to 8 hours of battery life and included charge case provides added on-the-go power



Water and dust resistant IP54 rating

2-year manufacturer's repair warranty

www.hphearingpro.com